# Duceralloy N

<b>D</b> Gebrauchsanweisung	<b>P</b> Instruções de uso	<b>CZ</b> Navod k použití
<b>GB</b> Instructions for use	<b>PL</b> Instrukcja stosowania	<b>TR</b> Kullanım Talimatı
<b>I</b> Istruzioni per l'uso	<b>SK</b> Návod na používanie	<b>H</b> Használati utasítás
<b>F</b> Mode d'emploi	<b>SLO</b> Navodila za uporabo	<b>UAE</b> إرشادات الاستخدام
<b>E</b> Instrucciones de uso	<b>RUS</b> Правила пользования	<b>ROK</b> 加工说明
<b>NL</b> Gebruiksaanwijzing	<b>GR</b> Οδηγίες Χρήσης	<b>CN</b> 加工处理 注意事项

**CE** 0124

**R<sup>x</sup>only**

**DeguDent**  
A Dentsply Company

<b>(D)</b>	Gebrauchsanweisung	1	<b>(SLO)</b>	Navodila za uporabo	55
<b>(GB)</b>	Instructions for use	7	<b>(RUS)</b>	Правила пользования	60
<b>(I)</b>	Istruzioni per l'uso	13	<b>(GR)</b>	Οδηγίες Χρήσης	66
<b>(F)</b>	Mode d'emploi	19	<b>(CZ)</b>	Navod k použití	72
<b>(E)</b>	Instrucciones de uso	25	<b>(TR)</b>	Kullanım Talimatı	77
<b>(NL)</b>	Gebruiksaanwijzing	31	<b>(H)</b>	Használati utasítás	82
<b>(P)</b>	Instruções de uso	37	<b>(UAE)</b>	إرشادات الاستخدام	87
<b>(PL)</b>	Instrukcja stosowania	43	<b>(ROK)</b>	加工说明	92
<b>(SK)</b>	Návod na používanie	50	<b>(CN)</b>	가공처리 주지사항	98

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## Venerable NiCr dental casting alloy

### Intended purpose:

Duceralloy N is a highly corrosion-resistant NiCr-based NPM deflagration alloy. It can be veneered with all dental ceramic masses covering the CTE range from 14.4  $\mu\text{m}/\text{m}\cdot\text{K}$  (25 – 600 °C) (e.g. StarLight Ceram, Ceramco 3, Ceramco iC, Duceram Plus, Duceram Kiss, Duceram love). Duceralloy N complies with the specifications of the standard DIN EN ISO 22674.

### Material data:

#### Composition:

Ni	61.0 %
Cr	25.8 %
Mo	11.0 %
Si	1.5 %
W	< 0.5 %
Al	< 0.4 %

Free of beryllium

Technical data:

Type:	Type 3 NPM deflagration alloy
Vickers hardness:	185 HV 10
0.2% expansion range:	340 MPa
Tensile strength:	550 MPa
Modulus of elasticity:	115 GPa
Fracture strain:	26 %
Density:	8.2 g/cm <sup>3</sup>
Solidus temperature:	1325 °C
Liquidus temperature:	1350 °C
CTE: (25 – 500 °C)	13.8 µm/m·K
(25 – 600 °C)	14.4 µm/m·K

**Contraindications:** Do not use in cases of confirmed hypersensitivity to one or more of the metals in the alloy.

**Warning:** Duceralloy N contains chromium and nickel. May induce allergic reactions.

This alloy should not be integrated or used in cases of known sensitization to one of the components.

**Precautions:** When handling dental alloys, please remember that the dust and vapours must not

be inhaled. Use suitable extraction equipment to protect persons from dust and vapours and wear face masks or other respiratory protection as well.

**Adverse effects:** Allergies may be induced to the metals in the alloy and electrochemical paraesthesias are also possible. Isolated cases are reported of systemic adverse effects of metals contained in the alloy.

**Interactions:** Avoid occlusal and approximal contact between different alloy types.

### Processing information:

**Modelling:** Cover stump with stump lacquer before modelling. The form of the crown should correspond to the reduced tooth. When modelling make sure not to create any sharp corners or edges.

They could subsequently cause cracks in the veneer ceramic. When modelling and working up the frameworks, make sure not to go below the minimum wall thicknesses of 0.3 mm for single crowns and 0.5 mm for pillar crowns.

**Pinning:** If bar casting is not preferred and single crown pinning is used instead, the casting channel should have a diameter of 3 – 4 mm. The following casting channel dimensions should be maintained if the bar casting method is used:

Casting channel from cone	3 – 4 mm diameter
---------------------------	-------------------

Casting channel crosswise (bar)	4 – 5 mm diameter
---------------------------------	-------------------

Connection to crowns	3 – 3.5 mm diameter and 5-8 mm in length
----------------------	--

**Investment:** Any phosphate-bound investment compound for high-fusing alloys that is free of graphite can be used (e.g. Deguvest SR). A sufficient muffle lining, if necessary with several layers, is required for work with metal muffle rings. Follow the instructions for the investment compound for the investment procedure.

**Wax removal / preheating:** At 300°C and 900°C depending on the size of the casting cuvette:

Casting cuvette size	1 x 20 min
	3 x 30 min
	6 x 45 min
	9 x 60 min

**Amount of metal required:** The amount of material required for the cast is calculated according to this rule of thumb:

"Wax weight of modelling multiplied by the alloy density (8.2 g/cm<sup>3</sup>) + approx. 10 g".

**Fusing method:** Duceralloy N can be melted using both the high-frequency method and the open-flame method.

**High-frequency method:** The cast is released 3-5 seconds after complete fusing.

**Induction method (vacuum casting method):** The cast is released 3-5 seconds after breaking of the oxide skin.

**Open flame casting:** To fuse Duceralloy N the burner must be used such that the heat necessary for fusing is evenly distributed. Use maximum burner setting. Follow instructions of burner manufacturer. The alloy forms an oxide skin when it melts. Continue application of heat until the flame pres-

sure causes visible movement of the metal under the oxide skin. The cast is released after 3-5 seconds.

**Important! Do not use melting powder.**

**Crucible materials:** Only ceramic crucibles may be used for casting Duceralloy N. The crucible should be preheated in the oven without metal before casting.

**Workup:** Carbide metal cutters are suitable for working up the casts. Diamond tools should not be used. Sharp tools are important in any case. Always work at high cutting speeds and low pressure levels.

**Oxide firing:** Oxidation is not necessary in principle.

**Ceramic veneering:** After workup, the surface of Duceralloy N is roughened by intensively blasting the surface with  $Al_2O_3$  (250  $\mu m$ ) at 3-4 bar.

Deflagration of the ceramic is done in accordance with the instructions for use provided for the ceramic compound used. Ideally suited are the ceramics StarLight Ceram, Duceram Plus, Duceram Kiss, Duceram love, Ceramco 3 and Ceramco iC. Comply with the recommendations of the ceramic manufacturer for NF alloys.

**Polishing:** To reduce oxidation of the unveneered alloy segments, prep polishing of the surfaces prior to ceramic firing is recommended. After firing, the alloy surface is worked over with a rubber polisher until a silky matte sheen is obtained, the polished to a high gloss with polishing or diamond paste. Polishing of the metal surface increases its corrosion resistance. This reduces the release of alloy components and thus further improves the biocompatibility of the alloy. Polishing agent residues can then be removed using a steam jet device or in an ultrasonic bath.

### **Soldering:**

Before firing: Degudent Solder U1W (1,120°C)

Flux: Oxynon.

### **Laser:**

Biosil laser welding wire

**Reuse of used material:** Sprue cones and channels or miscasts must not be reused since this results in oxide inclusions that compromise quality.

**For dental use only. Store under lock and key and out of reach of children.**

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DeguDent GmbH