## R-080110

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The R-080110 is a practical furnace especially designed to work with zirconium in the field of dental prosthesis. The unit features a control panel to program a sequence of different temperatures and times according to the user's needs.

## INSTALLATION

Carefully unpack the furnace. Inside the package you will find:

- 100 g balls.
- 1 crucible plate.
- 1 isolation block.
- Pliers.

Check out that the muffle that covers the inside of the furnace has not been damaged during transport.

### ▪ Installing the resistors:

1. Disconnect the power cord of the furnace from the power source.
2. Unscrew the 4 bolts in the top panel of the furnace and remove it (Fig. 1).
3. Insert the resistors into the slots (Fig. 2) together with the ceramic insulators (Fig. 3).
4. Adjust the distance to the lower part of the. Fix the height of each resistor by tightening the bolt of the clamp (Fig. 5)
5. Connect the terminals of the 4 resistors in series, as in Fig. 5.
6. Make sure that the resistors and the electric connectors are properly tightened.


Fig. 1

Fig. 2

Fig. 3

Fig. 4


Fig. 5

- Place the furnace on a horizontal surface, flat and firm, away from vibration sources. Make sure there is enough room for ventilation.
- Make sure the casters of the furnace are braked or locked during operation.
- The closest element should be at least 15 cm away from the unit.
- Connect the unit to a grounded 230 V, 50/60 Hz power source. It is also important that the power source is protected by a thermal magnetic circuit breaker of an appropriate voltage (consult the installation engineer).
- Before operating the furnace, it should be pre-heated in order to evaporate the humidity of the muffle. To do this, program the furnace to keep a temperature of 300 °C for an hour.

**BASIC PROGRAMMING**

The electronic controls of the R-080110 allow to program up to 50 steps. Each step consists of two values:

- C Temperature at which the furnace will start each step.
- t Time selected to reach the next temperature.


Fig 9

The parameters for the example graphic (Fig. 9) would be the following:

Step	Parameter	Value	Meaning
Step 1	C1	0 °C	Initial temperature
	t1	30 min.	Time selected to reach next temperature
Step 2	C2	400 °C	First set-up temperature
	t2	15 min.	Time selected to reach next temperature (maintenance)
Step 3	C3	400 °C	First set-up temperature (maintenance)
	t3	30 min.	Time selected to reach next temperature
Step 4	C4	1000 °C	Second set-up temperature
	t4	15 min.	Time selected to reach next temperature (maintenance)
Step 5	C5	1000 °C	Second set-up temperature (maintenance)
	t5	15 min.	Time selected to reach next temperature
Step 6	C6	600 °C	Third set-up temperature
	t6	-121	End

**Modifying the program:**

1. Switch on the circuit breaker (C) in the left side of the furnace.
2. Turn the knob (H) clockwise to switch on the control panel and modify the program.
3. Press the ◀ key. The upper display will show the [ 0 ] message.
4. Select the first set-up temperature by pressing the ◀▲▼ keys (the ◀ key allows to move from units to tens and to hundreds in order to modify their values more quickly). We recommend to select always 0 °C as initial temperature.
5. Press the ⏪ key. The upper display will show the [ 0 ] message
6. Select the first time by pressing the ▲▼ keys.
7. Press the ⏪ key. The upper display will show the [ 0 ] message.
8. Proceed in the same way for all the steps.
9. If you want to finish a working cycle, select the "-121" for the last time-parameter.
10. If you want to move directly to a given step, first press the ⏪ key; then select the step by pressing the ▲▼ keys, and access the parameters with the ◀ key.
11. If you want to leave the programming, press the ◀ key and the ⏪ key at the same time.

**Observations:**

- We recommend to use always the "0" value for C1.
- If you want to maintain the same temperature for a while, select that temperature in the next step.
- If you want to finish a working cycle, select the "-121" negative value for the last time-parameter.
- When you modify the value of a parameter, the value is automatically stored in the memory.
- The value of the parameters can be modified even when the program is under way.
- When selecting the time to reach a given temperature, select a time long enough to prevent the furnace from working at full power when it is not necessary.
- **Do not program a time lower than 15 minutes from 0 to 200 °C.**

**ADVANCED PROGRAMMING**

The electronic controls of the R-080110 allow to store different programs in the memory, one after the other, and to use each one whenever it is needed.

To do this, insert the values of the different programs from Step 2 on, and finish each program with the "-121" time value.

You can use the t01 value, with a negative sign, to select the step in which you want the furnace to start working.


Fig 10

If we want to store the programs of the example graphics (Fig. 10) and carry out the second program, the parameters should be the following:

Step	Parameter	Value	Meaning
Step 1	C1	0 °C	Initial temperature
	t1	-8	Step in which the furnace will start working
Step 2	C2	0 °C	Initial temperature
	t2	30 min.	Time selected to reach the next temperature
Step 3	C3	400 °C	First set-up temperature
	t3	15 min.	Time selected to reach the next temperature (maintenance)
Step 4	C4	400 °C	First set-up temperature (maintenance)
	t4	30 min.	Time selected to reach the next temperature
Step 5	C5	1000 °C	Second set-up temperature
	t5	15 min.	Time selected to reach the next temperature (maintenance)
Step 6	C6	1000 °C	Second set-up temperature (maintenance)
	t6	15 min.	Time selected to reach the next temperature
Step 7	C7	600 °C	Third set-up temperature
	t7	-121	End
Step 8	C8	0 °C	Initial temperature
	t8	15 min.	Time selected to reach the next temperature
Step 9	C9	500 °C	First set-up temperature
	t9	60 min.	Time selected to reach the next temperature
Step 10	C10	1500 °C	Second set-up temperature
	t10	20 min.	Time selected to reach the next temperature (maintenance)
Step 11	C11	1500 °C	Second set-up temperature (maintenance)
	t11	-121	End

## OPERATION

1. Put the material that you want to heat into the furnace. To do this, whenever the furnace is hot, wear gloves and use the pliers. TALLERES MESTRAITUA, S.L. strongly recommends to use **MESTRA** pliers.
2. Close the door and lock it with the key to prevent accidental opening.
3. Switch on the circuit breaker (C) in the left side of the furnace.
4. Turn the **LOCK** knob (H) clockwise to switch on the control panel and modify the program.
5. Program the furnace as explained in the previous section.
6. Pulse el botón verde (J) para conectar el circuito de potencia.
7. The lower display will show the "STOP" message, which means the furnace is not working
8. Keep the **▼/RUN** key pressed for a few seconds. The furnace will start up.
9. The upper display shows the real temperature inside the muffle.
10. The lower display shows the theoretical temperature the controls are trying to reach in each case.
11. Once the program is finished, the display will show the "STOP" message.
12. If you want to finish the program under way, keep the **▲/STOP** key pressed for a few seconds.
13. If you want to interrupt a program under way, keep the **▼/RUN** key pressed for a few seconds. The electronic controls will show the "Hold" message and interrupt the timer; the controls will try to maintain the theoretical temperature until the program is restarted again by keeping pressed the same key for a few seconds.
14. When you have finished the work, press the red button (I) to disconnect the power circuit.
15. **Do not open the door until the temperature of the muffle has dropped below 300 °C. A sudden change of temperature may damage the resistors.**
16. When the furnace is hot, wear gloves and use the pliers to remove the material from inside the muffle. TALLERES MESTRAITUA, S.L. recommends the use of **MESTRA** pliers.
17. To switch the control panel off, turn the **LOCK** knob (H) anticlockwise. If you are not going to use the furnace for a long time, switch off the circuit breaker (C) in the left side of the furnace.

The furnace has two thermal magnetic circuit breakers to protect it from electric overloads. One of them is placed in the left side of the furnace, and the other one inside it (Fig. 11). Once the problem has been solved, the circuit breakers have to be rearmed again.

**WARNING:** Connect the power circuit (green button) before you start the pre-heating cycle with the electronic controls (by keeping the **▼/RUN** key pressed). Otherwise, an overload could make the circuit breaker inside the furnace go off or even the lab electrical installation.

Fig. 11

## MAINTENANCE AND CLEANING

For a high reliability and a long working life of the unit, we strongly recommend the following:

- Clean the outside of the furnace periodically.
- When putting in or removing the materials from the furnace, take care not to damage the resistors or the muffle that covers the inside of the furnace.
- Prevent any liquid or solid substance dirtying or soiling the furnace muffle.
- Do not open the door when the temperature is above 300 °C. A sudden change of temperature may damage the resistors.
- Handle the resistors very carefully, since they are very fragile.
- When selecting the time to reach a given temperature, select a time long enough to prevent the furnace from working at full power when it is not necessary

**ERROR MESSAGES**

<i>Symptom</i>	<i>Reason</i>	<i>Solution</i>
The electronic controls is unblocked, but green pilot does not light up.	The circuit breaker inside the furnace has gone off (Fig. 11).	Check the furnace and rearm the circuit breaker.
The green light and the red light are turned on at the same time.	The circuit breaker inside the furnace has gone off (Fig. 11).	Check the furnace and rearm the circuit breaker.
The amperimeter does not show any strength, but the voltmeter shows its maximum value.	Resistor broken.	Find out the broken resistor and replace it.
The display shows the "OrAL" message.	Thermocouple broken	Replace the thermocouple.
The display shows the "HIAL" message.	The temperature in the furnace has gone beyond the maximum (1650 °C).	Let the furnace cool down and find out the reason for such a high temperature (the program may be wrongly designed).
The power source and the resistors are all right, but the temperature does not rise.	Electronic controls damaged.	Contact technical service.

**CAUTION**

- Before connecting the unit, make sure the power source is a grounded 230 V, 50/60 Hz. The power source should be properly protected by a thermal magnetic circuit breaker (consult the installation engineer).
- Do not allow children or unqualified personnel to handle the unit.
- Always keep the door closed while the furnace is working.
- Open the door with caution. Fumes coming out suddenly may cause burns.
- Do not load the furnace excessively, and do not insert objects that may fall when opening the door.
- Do not put flammable materials or liquids into the furnace or substances that may produce toxic fumes.
- Place the furnace on a flat, firm surface. Make sure it is away from vibration sources and flammable materials.
- Place the furnace where its sides are at a minimum distance of 15 cm from the wall or any other object.
- Wear gloves and use pliers whenever you insert or remove materials from the furnace. Do not touch the metal casing, since it may reach a high temperature.
- After heating the furnace, remove the materials with extreme caution. They may reach a very high temperature (up to 1650 °C) and cause accidents.
- Do not touch the resistors when the furnace is connected to the power source.
- Do not touch the resistors without making sure they are cold.
- Disconnect the power cord before carrying out any maintenance operation.
- Do not ever insert your fingers or any object through the ventilation grids

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