

PRODUCT OVERVIEW



Translucency - Added Value & Aesthetics - Shades Classification

priti[®]multidisc ZrO₂

p^{riti}® multidisc ZrO₂

Suitable for all conventional open milling machines with adapter for **blank diameter 98.5 mm**.

p^{riti}® multidisc ZrO₂ – High Translucent with a **translucency of 49%**, Sintering temperature 1450° C
Flexural strength > **650 MPa** - Highest Aesthetics - Released for bridges up to 3 units

HIGH TRANSLUCENT	SHADES	HEIGHTS IN MM
white	white	10, 12, 14, 16, 18, 20, 25
monochrome in 16 VITA® classical shades	A1, A2, A3, A3,5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	10, 12, 14, 16, 18, 20, 25
multicolor	A light (shade range: A1-A3), A dark (shade range: A3,5-A4), B light (shade range: B1-B2), B dark (shade range: B3- B4) C light (shade range: C1-C2), C dark (shade range: C3-C4) D light (shade range: D2-D4), p ^{riti} ® bleach (shade range: 0M1, 0M2, 0M3, 1M1)	14, 18, 20, 25

p^{riti}® multidisc ZrO₂ – Translucent with a **translucency of 40%**, Sintering temperature 1450° C
Flexural strength > **1,150 MPa** - High Aesthetics - Released for bridges up to 16 units

TRANSLUCENT	SHADES	HEIGHTS IN MM
white	white	10, 12, 14, 16, 18, 20, 25
monochrome in 16 VITA® classical shades	A1, A2, A3, A3,5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	10, 12, 14, 16, 18, 20, 25
multicolor	A light (shade range: A1-A3), A dark (shade range: A3,5-A4), B light (shade range: B1-B2), B dark (shade range: B3- B4) C light (shade range: C1-C2), C dark (shade range: C3-C4) D light (shade range: D2-D4)	14, 18, 20, 25

p^{riti}® multidisc ZrO₂ – Opaque with a **translucency of 35%**, Sintering temperature 1450° C
Flexural strength > **1,150 MPa** - Released for bridges up to 16 units

OPAQUE	SHADES	HEIGHTS IN MM
white	white	10, 12, 14, 16, 18, 20, 25
monochrome in 16 VITA® classical shades	A1, A2, A3, A3,5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4	10, 12, 14, 16, 18, 20, 25

All items can be ordered with or without a step (delivery time on request).

TRANSLUCENCY

**MADE IN
GERMANY**







	Opaque	Translucent	High Translucent
Multicolor		 priti®multidisc ZrO ₂ multicolor Translucent	 priti®multidisc ZrO ₂ multicolor High Translucent
Monochromatic	 priti®multidisc ZrO ₂ 16 VITA® shades A1 - D4 Opaque	 priti®multidisc ZrO ₂ 16 VITA® shades A1 - D4 Translucent	 priti®multidisc ZrO ₂ 16 VITA® shades A1 - D4 High Translucent
White	 priti®multidisc ZrO ₂ white Opaque	 priti®multidisc ZrO ₂ white Translucent	 priti®multidisc ZrO ₂ white High Translucent

After Added Value

After Aesthetics

ADDED VALUE & AESTHETICS


multicolor	
<p>High Translucent, cubic zirconia</p> <p>49% translucency flexural strength > 650 MPa bridges with up to 3 units</p> <p>multicolor</p>  <p>Variants multicolor</p> <p>1 bleach shade 4 light variants 3 dark variants</p>	<p>Translucent</p> <p>40% translucency flexural strength > 1,150 MPa suitable for a wide range of indications</p> <p>multicolor</p>  <p>Variants multicolor</p> <p>4 light variations 3 dark variants</p>

monochromatisch		
<p>High Translucent, cubic zirconia</p> <p>49% translucency flexural strength > 650 MPa bridges with up to 3 units</p> <p>16 VITA® shades and white</p>   <p>Variants monochromatic</p> <p>16 VITA® shades</p> <p>Variants white</p>	<p>Translucent</p> <p>40% translucency flexural strength > 1,150 MPa suitable for a wide range of indications</p> <p>16 VITA® shades and white</p>   <p>Variants monochromatic</p> <p>16 VITA® shades</p> <p>Variants white</p>	<p>Opaque</p> <p>35% translucency flexural strength > 1,150 MPa suitable for a wide range of indications</p> <p>16 VITA® shades and white</p>   <p>Variants monochromatic</p> <p>16 VITA® shades</p> <p>Variants white</p>

Why does pritidenta offer so many different translucencies, shades and heights?

Each lab or milling center has its own specific requirements regarding result, price, translucency and area of application: priti®multidisc ZrO₂ offers the right solution for every need - Good to know you have a choice when you need it!

VITA® and VITA® color designations are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.



pritidenta® GmbH


Meisenweg 37

D-70771 Leinfelden-Echterdingen
Germany

E-Mail info@pritidenta.com

Phone +49 (0) 711.3206.56.0

Fax +49 (0) 711.3206.56.99



www.pritidenta.com/pritimultidisc